# signum<sup>®</sup>

- DE** Produktinformation
- GB** Product information
- FR** Information Produit
- ES** Información de producto


**KULZER**  
MITSUI CHEMICALS GROUP

# signum<sup>®</sup>

- [IT](#) Informazioni sul prodotto
- [PT](#) Informações do produto
- [CZ](#) Informace o produktu
- [PL](#) Informacje o produkcie

**KULZER**  
MITSUI CHEMICALS GROUP

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
Signum composite, Signum composite flow und Signum matrix are light-curing hybrid composites for crown and bridge work.

Signum ceramis is a light-curing glass ceramic composite for framework-free, fixed restorations.

### **Areas of application for Signum composite/Signum composite flow/ Signum matrix**

#### **Suitable for veneering the following dental framework materials:**

- Metal alloys (precious metal, non-precious or titanium)
- Zirconia
- High-performance polymers (PEEK)

#### **Suitable for partial alteration of the colour and shape of the following dental materials:**

- PMMA
- Photopolymers

## **Suitable for the following prosthetic restorations:**

- Full and partial veneers for permanent framework-supported crowns and bridges
- Veneering removable combination prosthetics (telescopic and tapered crowns and attachment prosthetics)
- Veneering of superstructures/tertiary structures on implants
- Metal-free front tooth and posterior tooth temporary restorations
- Alteration of the colour and shape of Kulzer acrylic denture teeth
- Alteration of the colour and shape of acrylic-based frameworks (PMMA)
- Alteration of the colour and shape of Kulzer photopolymers that are approved for characterisation with Signum components (according to the photopolymer instructions for use)

## **Note**





The framework-supported full veneering using Signum composite/Signum composite flow/Signum matrix is approved up to the 2<sup>nd</sup> premolar.

As long as excessive discomfort is not anticipated in the molar region, Signum composite, Signum composite flow and Signum matrix can be used to achieve aesthetically demanding results, also for the occlusal veneering used in this region.

If unphysiological discomfort is anticipated (e.g. bruxism or parafunctions) in the occlusal surface in molars, the veneer is to be placed without contact or provided with stops.

## **Areas of application for Signum ceramis**

- Framework-free single crowns for anterior and posterior teeth
- Veneers, inlays, onlays, partial crowns
- Temporary anterior and posterior crowns
- Alteration of the colour and shape of Kulzer acrylic denture teeth

framework-supported

framework-supported

## Framework design and material thickness – Signum composite/Signum composite flow/Signum matrix

### – Material thickness

See the information and requirements in the instructions for use of the framework material indicated for veneering.

We recommend increasing the minimum thickness with softer alloys, mainly at the cervical to ensure stability. This is particularly important for telescope and tapered crowns.

- **The framework/Signum composite/composite flow/matrix transition sites** must be outside of the antagonist contact or contact point (partial veneering).
- For **full veneering with Signum composite/composite flow/matrix**, a minimum layer thickness of 1.2 mm occlusal (fissures) and 1.2 mm incisal is required. The cross section of the approximal connectors should be stable enough to resist distortion. If there is insufficient space occlusally, fabricate a stop with a large surface area so that all functional excursions are in the region of the stop.

### – Full/partial facing of implant superstructures

Do not place the framework/**Signum composite/composite flow/matrix** transition in occlusal screw connections in the functional region of the antagonists. The entire screw access channel is designed in the framework material.


If the framework structures are reduced, make sure there is functional support.

Before veneering, apply **Signum effect ET1** to the thin areas of the framework.


- **The anchor crown and pontics should only be 2/3 of the natural size of the teeth.**

### Note

We recommend **Kulzer** alloys. For more information, refer to the alloy table.

framework-free

framework-free

## Material thickness – Signum ceramis

The **framework-free Signum ceramis crown** requires a chamfer or shoulder preparation of 0.5 mm.


**Crown:** Ensure that the occlusal minimum thickness of 1.2 mm is achieved.

The lingual surface of the upper incisors must not be less than 1.2 mm.

**Inlay:** The minimum isthmus width of 2 mm must be achieved, occlusal minimum layer thickness: 1.2 mm

**Veneers:** Minimum layer thickness 0.5 mm (cervical), 0.7–1 mm (central-vestibular), 1.3–1.5 mm (incisal)

Wax up crown and bridge frameworks anatomically for single surface facings and affix retentions (beads 0.4mm) to metal frameworks.

Fully prepare the metal framework and polish; reduce the retention beads to the equator from 0.4 to 0.2mm. Sandblast the surface to be faced with  $Al_2O_3$  (110–125 $\mu$ m, 3 bar), clean with oil-free compressed air or clean tap water (do not use steam jets) and dry carefully. Apply Signum metal bond I with a disposable brush and allow to dry.

## Framework preparation

Surfaces of the framework that are not to be masked with Signum opaque F should be polished to a high-lustre before application.

Clean the framework and sandblast the surfaces to be veneered.

### Note

Sandblast the surface to be veneered with 110–125 $\mu$ m of  $Al_2O_3$  (disposable blasting media) (metal: 3 bar, PEEK: 2 bar, zirconium dioxide: 2 bar).

At a sandblasting angle of 45° the optimum distance between the sandblasting nozzle and framework surface is 2–4cm. Next, clean metal and PEEK surfaces with clean (oil-free) compressed air or clean tap water (no steam jets) and dry thoroughly. Clean zirconium dioxide surfaces using a steam cleaner, or clean tap water, and dry thoroughly with clean (oil-free) compressed air. Avoid surface contamination.

## Framework conditioning

Apply Signum composite/composite flow/matrix/Signum opaque F for veneers on frameworks in conjunction with verified bonding systems Signum metal bond or signum universal bond (Refer to Signum metal bond and/or Signum universal bond instructions for use).

Apply a layer of Signum metal bond II and polymerise.

Fill the pontic with Signum composite effect ET 1 and polymerise (refer to the table on Page 17).

Apply the first and second layer of Signum opaque F and polymerise each layer separately.

## Signum metal bond or Signum universal bond

Signum metal bond I or Signum universal bond I is now placed drop by drop, per veneer unit, in a serving bowl and applied to the dry, dust- and grease-free surface with a disposable brush. Allow to air dry.

Shake Signum metal bond II several times in order to achieve an optimal distribution of the color pigments on the surface during application. Then Signum universal bond II is filled into a serving bowl and applied with another disposable brush – cure in a polymerisation device (refer to the table on page 17 for polymerisation times).

## Signum opaque F

### Note

For pontics, apply a layer of Signum opaque F and polymerise it prior to application of Signum composite effect ET1. Fill Signum composite effect ET1 to the level of the adjacent crowns (max. 5mm). Apply Signum opaque F to the polymerised Signum composite effect ET1 as described.

Signum opaque F is placed on the mixing block and stirred thoroughly. Apply the opaque with a suitable, disposable brush in several, uniform, semi-overlapping layers until the scaffolding frame is completely covered. Each layer must be polymerised individually. To obtain an optimum adhesive bond, it is necessary to apply the opaque without contaminating the intermediate layer (dispersive layer) caused by touching it or by long pauses. Next, apply the components of the veneering system from Kulzer.

Apply a layer of Signum composite margin/Signum composite flow cervical to the cervical crown margin and polymerise.

Apply Signum dentine and polymerise.

Apply and polymerise Signum composite effect or Signum cre-active for customising the shade (cover Signum cre-active masses with Signum composite).

Apply Signum composite enamel in the upper third of the build-up and polymerise. Prepare the veneering using the Signum tool kit and polish to a high lustre.

## Signum composite/composite flow

Before building up further, separate all areas on the stone model that may come into contact with Signum components. The composite veneer is designed in layers in a 3-layer technique (margin-dentine-incisal-layer) or, depending on the desired aesthetic result, with additional individualised dimensions (Signum matrix), corresponding to their respective application areas. Apply the individual components only on layers already polymerised. The composite only achieves its physical and colour properties by the polymerisation of each individual layer.

Apply Signum pastes in portions.

Condense the Signum paste dispensed from the cartridge with a sculpting instrument, apply to the polymerised Signum layers and press firmly. Contour with a sculpting instrument or brush. Condensing the Signum paste prevents air bubbles becoming trapped between the layers.

### Note

If the Signum composite enamel layer is too thick or the dentine core is built up too much, it makes the shade of the facing appear grey. The optimum shade is reproduced with a layer thickness of 0.5 mm–2 mm. The PALA cre-active gum mass is available for supplementing gingival portions.

Finished veneering with components from the Signum essentia Set.

### Signum matrix

Low viscosity Signum matrix can be used instead of the high viscosity paste for the incisal area. These composites can be pregelled using der HiLite pre 2 to facilitate application prior to polymerisation (refer to the table on page 17 for polymerisation times).

The finished margin/dentine build-up

The value build-up


Apply the value composites (VL1 – VL4), which are highly fluorescent, in order to adjust the brightness.

The correct brightness value of thinner facings with dark frameworks showing through is also very easily attained using value composites.

Mamelon dentine and secondary dentine

Mamelons can be applied with a brush or sculpting instrument. Secondary dentine can be used if required for the cervical area, palatal surfaces or imitating abraded areas. The mamelon and secondary dentines MD1, MD2, MD3 and SD1, SD2 are also suitable for use as opaque dentine due to their lower transparency.

Briefly stir the OS and OT composites with an instrument or brush to obtain a more fluid consistency for a short period.

This can be an advantage when the Signum matrix components need to be more flowable on application.

Opal transpa build-up

The opal transpa composites (OT1–OT10) are applied to the incisal area to reproduce the opalescence and transparency of a natural tooth.

The smaller the number, the greater the transparency effect of the opal transpa composite. OT1 is the composite with the highest transparency and lowest opalescence, while OT10 has the lowest transparency and highest opalescence.

In addition to adjusting the transparency and opalescence, the coloured opal transpa composites OTA, OTB, OTG, OTY, and OT ice can be used for characterising the shade.

Opal incisal build-up

Finally the anatomical shape is completed with an appropriate opal incisal (OS1–OS4) (also refer to the table on Page 16).

Prepared and polished  
Signum matrix restoration**Note**

The OS and OT, MD and SD composites can be pregelled using the HiLite pre 2.

Do not exceed a maximum layer thickness of 1 mm when pregelling (refer to the table on Page 17 for polymerisation times).

### Signum ceramis – Inlay/Onlay

Fabricate a sectional model using neutrally coloured high-strength or standard dental stone. Block out undercuts with a neutral colour of wax. First, separate the cavity and mark the preparation margin with a pencil.

Preparing the cavity

Apply a thin layer of Signum ceramis effect, e.g., ET1, to the cavity walls extending to the preparation margin to improve the shade adaptation with inlays/onlays. Apply translucent Signum ceramis effect composites, e.g., ET4, to the cavity floor to enhance the depth effect.

Build-up with Signum ceramis dentine, and Signum ceramis effect

Apply Signum ceramis dentine with a sculpting instrument into the cavity and contour it. Signum ceramis, Signum cre-active und Signum matrix components can be applied to polymerised Signum ceramis dentine for customised characterisation of the inlays/onlays.

Characterisation using Signum cre-active composites

Contouring the anatomical shape

Apply translucent Signum ceramis effect composites to the relevant areas and polymerise. The translucent Signum ceramis effect composites (ET1) allow the natural tooth shade to be transmitted to the Signum ceramis restoration at the contact points with the natural teeth (chameleon effect).

After surface treatment and polishing

After trimming and polishing (refer to trimming/polishing) with the Signum tool kit, carefully roughen the restoration base with a coarse diamond bur or 50µm Al<sub>2</sub>O<sub>3</sub> at 2 bar pressure (using disposable blasting media). Adhesive luting systems suited for this indication should be used to fix the restoration.

## Signum ceramis – single crowns

Fabricate a sectional model using neutrally coloured high-strength or standard dental stone (Octa Rock). Block out undercut areas as required. Cover the die with a thin layer of wax as a spacer. First, separate the prepared die and mark the preparation margin with a pencil.

Preparing the model/  
Separating the die

Signum cre-active or Signum matrix can be used for customised characterisation both in the cervical area and on the occlusal surface. These are polymerised before applying other Signum ceramis composites (refer to the table on Page 17 for polymerisation times). Apply Signum ceramis dentine and contour to the anatomical shape. For crowns in the posterior region, the buccal and lingual tooth cusps and marginal ridges are contoured first.

Characterisation

Apply Signum ceramis enamel to the upper third of the facing on the dentine and contour to the full anatomical shape with a sculpting instrument or a brush. Signum ceramis effect or Signum matrix can be used as well as Signum ceramis enamel for customising the shade of the restoration in the incisal and occlusal areas. Apply the relevant composites to the required areas and polymerised. Highly translucent composites (ET1, OT1) at the contact points to the natural tooth allow the tooth shade to be transmitted to the crown (chameleon effect).

Shade reproduction

After final polymerisation, remove the crown from the stone die and prepare it. To ensure optimal preparation of Signum ceramis, always use tungsten carbide cutters. (Refer to the section on trimming/polishing on Page 14 for details).

Finished restoration

After trimming and polishing (refer to trimming/polishing) with the Signum tool kit, carefully roughen the restoration base with a coarse diamond bur or 50 µm Al<sub>2</sub>O<sub>3</sub> at 2 bar pressure (using disposable blasting media). Adhesive luting systems suited for this indication should be used to fix the restoration.


## Signum connector

Light-curing Signum connector is used for bonding light-curing Signum composites to Kulzer denture base materials, Kulzer acrylic denture teeth or for characterising approved Kulzer photopolymers. Roughen the surfaces to be conditioned with a coarse-grained diamond or sandblast with  $Al_2O_3$  (at least  $50\ \mu m$  with a pressure of 2 bar). Blow off coarse dust with oil-free air. Fine dust particles that remain on the surface can affect the bonding and should be removed carefully with a clean brush with short, stiff bristles. At this stage, the oil and grease-free surface should not be steam-cleaned or touched. Put some Signum connector in a clean dish, apply a thin layer onto the clean surface with a disposable brush, allow to react for 2–3 mins and then polymerise (90sec.). Do not touch the surface after polymerisation. If the dispersion layer has been destroyed through wiping or sanding, the product application described here must be repeated. If the denture bases change color or shape, characterise them with Pala cre-active components. Artificial acrylic denture teeth are customised with Signum cre-active and/or Signum veneering material. Apply and polymerise the materials according to the recommended procedure (see the relevant instructions for use).

### Note

When bonding cold-curing resins used as denture base materials to Signum components, apply Signum liquid. Build-up occurs as previously described for Signum connector. If the substance is not allowed to take effect for at least 2 minutes, it will not create an adequate adhesive bond.

**Polymerisation:** 90sec. in Heraflash, HiLite power or HiLite power 3D

## Signum cre-active

The Signum cre-active composites are used for individual characterisation of light-curing Signum composites, Kulzer acrylic denture teeth and PMMA, as well as suitable Kulzer photopolymers. The ready-to-use Signum cre-active can be applied directly onto the polymerised Signum material. Signum cre-active components must not be mixed with Signum paste components to guarantee adequate polymerisation.

If Signum cre-active colorfluids have been applied, they must be covered with veneering material (enamel, opal incisal, opal transpa, etc.). If used on acrylic denture teeth, acrylic dentures or suitable Kulzer photopolymers, Signum connector should be used as the bonder.

### Note

All Signum cre-active components can be mixed with one another!  
PALA cre-active gingiva (flow) composites should not be used for veneering.

## Trimming/Polishing

The veneering can be prepared immediately following final polymerisation (refer to the table on Page 17 for polymerisation times). To ensure optimal preparation of Signum restorations always use fine cross-cut tungsten carbide cutters. Prepare buccal surfaces with Magnum and Piccolo cutters at approx. 15,000rpm. Smooth the framework/ veneering junctions with a Silico polisher (5–10,000rpm). Prepolish the facing with the Prepol (3,000–5,000rpm) without affecting the contour. This produces a silky sheen. Use the Mepol I with Signum HP paste or HP diamond and then the Hipol to polish to a high lustre. We recommend using the Mepol II for polishing Signum matrix and Signum ceramis to produce a durable high lustre. After cleaning the restoration with a steam cleaner, polish the facing again with a clean Hipol.

## Signum liquid

To make it easier to model or adapt light-curing veneering materials, the modelling tool or the brush are coated sparingly with Signum liquid.

To restore a damaged dispersion layer, the surface to be conditioned is roughened with a coarse grain diamond. Blow off coarse dust with oil-free air. Fine dust particles that remain on the surface can affect the bonding and should be removed carefully with a clean brush with short, stiff bristles. At this stage, the oil and grease-free surface should not be steam-cleaned or touched. Coat the point of attachment with Signum liquid applied in a thin and broad layer using a disposable brush. Make sure to moisten the area thoroughly. Allow to take effect for 30 seconds, do not polymerise. Apply the desired Signum components as usual.

## Polymerisation

Refer to the table on Page 17.

### Note

To bind denture base materials to light-curing veneering materials, prepare the area as described above. To create an adequate adhesive bond, allow the signum liquid to take effect for at least 2 minutes and then polymerise (90seconds in Heraflash, HiLite power or HiLite power 3D). Do not touch the surface after polymerisation. If the substance is not allowed to take effect for at least 2 minutes, it will not create an adequate adhesive bond.

## Signum metal bond/Signum universal bond

Condition the framework surface for facing as described in the relevant instructions for use. Signum metal bond I/Signum universal bond I is now placed drop by drop, per veneer unit, in a serving bowl and applied to the dry, dust- and grease-free surface with a disposable brush. Allow to air dry.

Shake Signum metal bond II several times in order to achieve an optimal distribution of the color pigments on the surface during application. Then Signum metal bond II/Signum universal bond II is filled into a serving bowl and applied with another disposable brush – cure in a polymerisation device (90s in Heraflash, HiLite power or HiLite power 3D). (Refer to the table on Page 17 for polymerisation times).

In order to achieve an optimal bond it is necessary to apply the opaque without contamination of the intermediate layer (dispersion layer) through contact or a longer break and then to apply the components of the veneering system or the denture base materials (hot or auto-polymerisate) from Kulzer.

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For aesthetic adjustments, a combination of all incisal build-ups with ET1 is recommended.

MD1; MD2; MD3; SD1; SD2; OT1–OT10; OTY; OTB; OTA; OTG and OT ice are used individually with all shades without a specific shade combination.

Signum opaque F

Signum composite flow

Signum matrix

Signum ceramis

Signum cre-active


Signum composite

## Polymerisation devices/times

Materials	Max. polymerisable layer thickness (mm)	HiLite power 3D, HiLite power, Heraflash (sec)
Signum metal bond II/Signum universal bond II	–	90
Signum connector	–	90
Signum opaque F	0.08	90
Signum composite effect ET1 For filling pontics	5.0	180
Signum composite margin	2.0	90
Signum composite dentine	2.0	90
Signum composite enamel	2.0	90
Signum composite effect	1.0	90
Signum composite flow dentine	2.0	90
Signum composite flow cervical	2.0	90
Signum matrix opal incisal OS	1.0	90
Signum matrix opal transpa OT	1.0	90
Signum matrix value VL	1.0	90
Signum matrix mamelon dentine MD	1.0	90
Signum matrix secondary dentine SD	1.0	90
Signum cre-active	0.2	90
Signum ceramis margin	2.0	90
Signum ceramis dentine	2.0	90
Signum ceramis enamel	2.0	90
Signum ceramis effect	2.0	90
Final polymerisation	alle 2.0	180*
<b>Materials:</b> margin, dentine, enamel, effect, value, opal incisal, opal transpa, mamelon dentine, cre-active, flow dentine, flow cervical. <b>Gelling units:</b> HiLite power 3D, HiLite power, Heraflash (6 sec.), HiLite pre 2 (10 sec.) Connect an intermediate polymerisation (90 sec.) after building up several layers up to the specified max. layer thickness.		

\* For injection and pressing technique: 2x 180 sec.

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Kulzer GmbH