MAINTENANCE, CONTROL AND CHECK

The products do not require special maintenance. Re-assemble disassembled products / individual components. You must always carry out a visual inspection for contamination, damage, wear and deformation before and after the individual work steps. Damaged or corroded products must no longer be used. If the product / individual component is not visibly clean, the entire processing must be repeated or the product / individual component must be disposed of properly.

PACKAGING

The product is to be packed in suitable and standardized transparent sterile packaging (sterilization bags) and sealed. Observe the instructions of the sterilization bag and sealing machine manufacturers as well as the current, normative requirements. Products / individual parts that have not been sterilized in bags must be used immediately.

STERILIZATION

Only tested steam vacuum autoclaves may be used. Make sure that during the sterilization of several different products / individual components, the autoclave is not overfilled and that the products / individual components do not contact each other.

The following sterilization cycles may be completed:

Steam sterilization, 134 °C, hold time 5 minutes or steam sterilization, 121 °C, hold time 15 minutes

To dry the products / individual components, the drying cycle of the autoclave should be set. Use the autoclave manufacturer's instructions for use.

STORAGE

To maintain the sterility, the products must be stored in standardized sterilization bags at a dry, clean place until they are used. If the sterile packaging is damaged, the products must be processed once again before their use.

ADDITIONAL INFORMATION

Mark the sterilized products according to the legal and national regulations. The recommended storage duration for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.

The user must ensure within the scope of their quality management system that specified processing cycles that might apply (see Limitations in the processing) are not exceeded.

The manufacturer and the competent authority of your member state shall be immediately notified of all serious incidents occurring in connection with the product.

Recommendations and regulations for disposal of the products vary by region. Ask your competent disposal company for the current regulations in your region.

CONTACT TO THE MANUFACTURER

Phone: (+49) 203 / 99269 -0

Fax: (+49) 203 / 299283

Email: info@hagerwerken.de

The instructions specified above are validated by the medical device manufacturer as suitable for the preparation of a medical device for its reuse. The processor is responsible for ensuring that the completed processing with the equipment, materials and personnel used in the preparation facility achieves the desired result. To this end, verification and/or validation and routine monitoring of the procedure are required.

GB

Instructions for Use

Inserting the electrodes into the handpiece

- Turn the head of the handpiece one turn to the left.
- Insert the shaft of the electrode into the handpiece as far as it will go.
- Turn the head of the handpiece to the right.

Now the electrode should fit tightly in the handpiece!

- Check the tightness of the electrode by gently pulling it. It must be impossible to turn the electrode or to pull it out of the shaft.
- Make sure that the non-coated part of the electrode (metal part) is inserted completely in the handpiece. No metal parts of the shaft must be visible outside the head of the handpiece. Otherwise, you pose a risk of injury to you and your patients.

IMPORTANT: Only the approved original electrodes of Hager & Werken must be used. These electrodes are adapted to the handpieces and devices.

Plugging the handpiece cables into the units

The handpieces are equipped with yellow or blue cables and matching plugs. Insert the plug into the corresponding socket with the same color cable and cable plug. Ensure that the plug is inserted quite deep - all the way to the stop.

The following applies to the hf1 Surg bipolar only

Always put both handpieces into the color-coded sockets and corresponding receptacles. This is the only way to use the device. By removing one handpiece from the receptacle, the second handpiece becomes functionless.

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GB

Processing instructions according to EN ISO 17664:2021

Product: HF handpieces
REF 452 423 (yellow), 452 425 (blue)

Manufacturer: HAGER & WERKEN GmbH & Co. KG, Ackerstraße 1, 47269 Duisburg, Germany

Warnings: Attention! The product must be cleaned, disinfected and sterilized according to these processing instructions before its first use as well as after every subsequent use. Disinfection alone is not enough. Thorough cleaning and disinfection are decisive prerequisites for effective sterilization. Processing should be started as soon as possible but not later than 2 hours after use. During processing, the product should not be unnecessarily exposed to wetness or humidity. Please also observe the legal regulations applicable in your country as well as the hygiene regulations of the medical practice and/or the hospital.

Limitations in processing: Frequent processing has only minor effects on this product. The end of the product's useful life is mainly determined by wear and damage caused by use. In case of doubt, the product should always be sorted out and replaced early. The user alone is responsible for making the decision regarding repeated use of the product. If the product is used too frequently, the manufacturer does not accept any warranty for its function, performance and safety.

Personal protective equipment: For reasons of occupational safety and to minimize infections / cross infections, suitable personal protective equipment (protective clothes, protective gloves, protective goggles and surgical masks) must be worn throughout the entire processing procedure.

INSTRUCTIONS:

INITIAL TREATMENT AT THE PLACE OF USE

Remove surface contamination using a disposable cloth / paper cloth when you are still at the place of use. Rinse the products with water (at least drinking water quality) 2 hours after the application at the latest. Avoid drying of residues or contamination of any kind on the product. When using the product for the first use, do not use any agents containing aldehyde or alcohol as they may cause protein fixation.

PREPARATION BEFORE CLEANING

Disassemble the product into its individual components. Perform a visual inspection for damage and wear.

CLEANING: MANUAL

Manual cleaning should be completed in a water bath (at least drinking water quality) with the cleaning agent specified, below using a brush below the water surface to achieve both sufficient cleaning of the products / individual components avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the degree of contamination of the product / individual component; it should, however, not be less than 1 minute. Afterwards, rinse the products / individual components with water (at least drinking water quality) for at least 20 seconds.

Description of the recommended cleaning process:

Pre-cleaning: Remove surface contamination using a disposable cloth.
Cleaning: Insert the instrument into the strainer bowl. Hang the strainer bowl into the cleaning bath with cleaning solution, remove contamination with cleaning brushes.
Cleaning tools: Flexbrush brush, REF 605 254
Rinsing: Rinse with water (at least drinking water quality) for 20 sec.
Drying: Dry at room temperature.

Manual cleaning should not exceed a temperature of 45 °C.

Ensure that visual contamination has been completely removed. If after the cleaning process, contamination is still visible, the cleaning has to be repeated.

DISINFECTION: MANUAL

For manual disinfection, please only use suitable disinfectants with tested effectiveness (CE mark, VAH-/DGHM-listed). Put the products / individual components into corresponding disinfection baths according to the disinfectant manufacturer's specifications. Ensure that the products are sufficiently covered and do not contact each other. The instructions for use of the disinfectant manufacturer must be strictly observed. In particular, observe concentrations and contact times.

Application with the following agents is recommended: - Dürr Dental ID 213 instrument disinfection.

Description of the recommended disinfection process:

Completion of proper cleaning preparation and manual cleaning.

Disinfection: Insert the part of the instrument that comes into contact with a patient into the strainer bowl. Hang the strainer bowl into the disinfection bath with disinfection solution.
Disinfectant: ID 213, company Dürr Dental, concentration: 2 %
Contact time: 5 min
Rinsing: Rinse with water (at least drinking water quality) for 20 sec.
Drying: Dry at room temperature.

Make sure that the disinfectant is compatible with any products and cleaning agents used. The pH of the disinfectant should be between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic / halogen hydrocarbons. Please observe thermal restrictions that might apply considering the disinfectant manufacturer's specifications. Afterwards, rinse the products / individual components with water (at least drinking water quality) for at least 20 seconds. Before further processing, the product must be free from any residues and dry.

DRYING

If the products are processed manually, the products / individual components may be dried by blowing them off with filtered, oil-free compressed air according to DIN ISO 8573-1 (medical cleanliness class) or at room temperature.