hf Surg[®]

HF-Chirurgiegerät

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HF-surgery device

User Manual

HF-chirurgisch apparaat

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Dispositif chirurgical HF

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Manual de instrucciones

Inhalt

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Index

Page 12	Advices and warranty / Scope of delivery
Page 13	Comments to content / Radio frequency (HF) – mode of action
Page 14 / 15	Technical data / Protection and safety regulations
Page 16	Preparation / Start the device / Lockout after use
Page 17	Labelling / Device description
Page 18	Performance graph
Page 19	Select the right HF-Mode / Practiac exercise on a beef mode / Example of clinical applications
Page 20 / 21	Reprocessing / Use Manual of hf handpieces
Page 52	Environment-directive

Please carefully read this user manual!

Read completely this user manual and get familiar with the use and functions of the unit and all accessories, before you start using the device clinically.

Unless you do not follow the directives as instructed, the following problems may occur:

- Serious injuries to the patient
- Serious injuries to the operator or to the service personal
- Damage of malfunction of the unit or of the accessories

Intended Purpose

The high frequency device hf-Surg is used for surgical applications in soft tissue in dentistry and is suitable for the following applications: oral surgery such as cutting, removal of soft tissue. No essential performance features according to EN 6061-1 are assigned to the device.

Field of application

The area of application is professional dentistry facilities.

Modifications

The manufacturer has the right to modify the appearance and technical data because of new product developments.

The marks: “WARNING“, “ATTENTION“ and “REMARK“ contain important hints.

Responsibility of manufacturer

Warranty and liability by Hager & Werken GmbH & Co KG is given, if:

- installation and start of operation is done by own personal or by personal authorized by the manufacturer
- installation and safety measures comply with national norms and regulations
- the unit is used in accordance with the user manual
- no manipulations to the unit or to the accessories, if manufacturer does not agree.

Warranty

The hf Surg unit has a legal warranty of 12 months.

Scope of Delivery hf Surg - REF 452 400

HF set of instruments:

Cutting electrode No. 1	REF 452 403
Cutting electrode No. 2	REF 452 404
Cutting electrode No. 15	REF 452 407
Cutting electrode No. 13	REF 452 411
HF hand piece Yellow 1,50 m	REF 452 423
Neutral electrode	REF 452 421
Foot pedal	
User manual	
Medical product journal	

Explanations to the user manual

Important directives, especially for technical safety and security, are mentioned:

ATTENTION

This information advises to special service procedures or caution measures, which must be considered to avoid damage to the unit.

REMARK

This is general and special information to clarify important and helpful instructions.

In the use of radio frequency devices, the emission of radiation is natural and cannot be avoided.

Interaction of radio frequency (HF)


If radio frequency is guided through very thin metal electrodes, a very large electro-magnetic power density is created in the tissue layer. The water is abruptly heated up in these cells which leads to cooking and rupture, respectively. High-frequency current is used to protect the patient from electric shock.

THERMAL INTERACTION (MONOPOLAR)

CUTTING (CUT) / COAGULATION (COAG)

At monopolar cutting the radio frequency current is led from the device via a work electrode and a large neutral electrode back to the device. The current density at the treatment point is very high, but it is very low at the neutral electrode.

The electric current density leads to a fast and strong heating in tissue, which can be used for cutting and coagulation.

This user manual contains information, which are copyright protected and are not allowed to be partly or completely copied or stored on microfilm or electronic items and should not be distributed without the written permission by Hager & Werken GmbH & Co KG.

If errors occur or if contents of this user manual are not clearly formulated, please inform Hager & Werken GmbH & Co KG.

Copyright © 2023 - Hager & Werken GmbH & Co KG

Technical Data

HF	frequency + power:	2,2 MHz, max. 50 W, load resistance 1000 Ohm
	method of operation	monopolar
	wave form	permanently / pulsed with ca. 80 %
Ambient temperature		+5 °C – +40 °C
Storage temperature		-20 °C – +70 °C
Humidity		Operating humidity: < 85%, non-condensing Storage humidity: < 90%, non-condensing
Air pressure		during storage : 500 hPa - 1080 hPa during operation : 700 hPa - 1080 hPa
Power supply		230 V AC, 50/60 Hz
Power consumption		max. 100 VA
Medical device classification		class IIb
Applicators	Hand piece	autoclavable
Safety		According to DIN EN ISO 6060-1
Electromagnetic compatibility		According to DIN EN ISO 6060-1-2
Safety hf surgery equipment		According to DIN EN ISO 60601-2-2

Cables

HF handpiece cable, yellow	Length 1.50m
HF neutral electrode cable	Length 1.50m
Foot switch cable	Length 2.50m
Power cable	Length 2.50m

Protection and safety provisions in the dental office for the use of the hf Surg device

The hf Surg device may only be put into operation after the instructions have been given by the operator and in compliance with the prescriptions and safety provisions.

WARNING

The device must only be connected to a power supply with an earth wire in order to avoid the risk of an electric shock. **Check the cables, hand piece and electrodes as well as the foot pedal on visible damage before starting up the hf Surg device. Instruments with brittle or faulty insulation must not be used as they pose a danger of injury.** When operating the hf Surg device, unpredictable malfunctions might occur that could cause unwanted output power increase.

1. During operation, a minimum distance of c. 20 cm from any walls must be kept.
2. The NEUTRAL ELECTRODE should be reliably stuck to a correspondingly prepared suitable area of the PATIENT body with its entire surface, as defined by the PRODUCER.3. The PATIENT should not get in touch with any metal parts that are grounded or have a significant earth capacitance (e.g. operating table supports etc.).
3. The PATIENT should not get in touch with any metal parts that are grounded or have a significant earth capacitance (e.g. operating table supports etc.).
4. Skin-to-skin contacts (e.g. between the arms and the body of the PATIENT) should be avoided, e.g. by inserting dry gauze. The cable leading to the hf Surg device should neither touch the patient nor any other lines. Instruments which are temporarily not in use during operation, a minimum distance of c. 20 cm from any walls must be kept.
2. The NEUTRAL ELECTRODE should be reliably stuck to a correspondingly prepared suitable area of the PATIENT body with its entire surface, as defined by the PRODUCER.3. The PATIENT should not get in touch with any metal parts that are grounded or have a significant earth capacitance (e.g. operating table supports etc.).
3. The PATIENT should not get in touch with any metal parts that are grounded or have a significant earth capacitance (e.g. operating table supports etc.).
4. Skin-to-skin contacts (e.g. between the arms and the body of the PATIENT) should be avoided, e.g. by inserting dry gauze. The cable leading to the hf Surg device should neither touch the patient nor any other lines. Instruments which are temporarily not in use during the treatment have to be kept away from the patient, e.g. on the instrument table.
5. The output power should be set as low as possible for the corresponding purpose.
6. An obviously low output value or functional failure of the hf Surg device in usual operation can be caused by insufficient adherence of the NEUTRAL ELECTRODE or insufficient contact in its connections. In such case, the adherence of the NEUTRAL ELECTRODE and its connections should be checked before setting a higher output power.
7. The use of ignitable anesthetic agents or combustible gases like nitrous oxide (N₂O) and oxygen should be avoided if a surgical intervention is executed in the head area, unless these substances are aspirated. If possible, non-ignitable ingredients should be used for cleaning and disinfection. Ignitable ingredients used as cleaning and disinfection agents or as solvents for adhesives should have evaporated before the use of the hf Surgery.
8. For patients with pacemakers or other active implants, there is a potential DANGER of disturbance of the pacemaker function or damage to the pacemaker. In case of doubt, an expert should be consulted.
9. The accessories must have a minimum accessory reference voltage of 500 V. Only use the original accessories contained in the delivery and offered by the producer to achieve maximum safety for the patient and the caregiver. The characteristics of the applied parts and conducts are adapted to the output power and output voltage of the device, so that a safe operation is ensured for all operation modes and settings.
10. The device must be disconnected from the power supply during cleaning.
11. Service and maintenance tasks may only be executed by authorized specialized personnel.
12. The radiation emitted by the hf Surg device during operation can interfere with the functionality of other electric devices. Computers, laptops and mobile phones should be kept away from the hf Surg device. Data on computers and laptops should be saved beforehand.
13. If any safety provisions and operating instructions contained in this manual are violated, any warranty and liability of the manufacturer is cancelled.
14. In operating rooms, the device may only be used with pedal switches with AP labeling.

Operator regulation

The device unit is classified as medical device unit class IIb (Europe).

Technical controls:

The user is committed to perform on a regular basis technical controls after the following specifications:

Period: every 24 months, starting with date of delivery and after each repair.

Covering:

Visual check of the unit and accessories - Check according to IEC 62353

- Protective earth resistance
- Alternate leakage current
- Alternate patient leakage current

Function check

- Main switch
- Switch CUT - CUT/COG with LEDs
- Uniformity of performance throughout the range of adjustment*

Measurement of radio frequency (HF) output power at a load of 1 k Ω :

- Output CUT (50W)
- Output COAG (45W)

All results of measurements must be documented according to DIN EN ISO 62353 concerning the first measured values. If defects occur during the controls, the user is responsible to initiate repair.

* Note for the user: the sound that can be heard during the operation of the HF changes its tone color with the operating mode (CUT or CUT/COAG), and the volume with the power setting. It is recommended that the user should check the regularity of the power output once in a while by turning the power knob from minimum to maximum with activated output and simultaneously checks if the signal has any dropouts. When changing the operating mode, the tone color must change a little bit. If dropouts or no difference between CUT and CUT/COAG can be heard, the device must be sent to the maintenance service.

Preparation for commissioning

- Before the hf Surg device is put into operation, it should be held at room temperature for a longer time (min. 30 minutes) to avoid the formation of condensation.

Commissioning

1. Insert the device plug on the back into the corresponding socket (n° 2) and the power plug into the power socket and connect the pedal switch (n° 3). Push the main switch (n° 1) of the device to the right, and one of the small yellow lamps on the front (n° 5 or 6) will turn on. Now, the device is ready for operation.

The small lamps in the yellow labeling field show the selected wave type:

yellow light above (n° 5) – cutting shaft

yellow light below [n° 6) – cutting/coagulation shaft B

By pushing the switch (n° 4), the wave type can be selected.

2. On the front part, connect the handpiece to the yellow socket
3. Connect the neutral electrode to the socket provided for that purpose (n° 3).
4. Put the desired electrode on the handpiece.
5. Check the intensity regulator (n° 1) and adjust it if required.
6. Activate the pedal switch, the pure tone typical of the set wave type can be heard.

Important notice:

Always activate the electrode by pushing the pedal switch before touching the tissue when you want to cut or cut/coagulate.

Lockout after use

1. Place the hand piece into the instrument holder
2. Remove the electrode and reprocess it accordingly
3. Switch off the device to position 0, using main switch 1 (back).

Label

At the back side of the unit:


Explanations

- ① Manufacturer
- ② On / Off
- ③ Mains fuse
- ④ Foot control
- ⑤ CE mark
- ⑥ Alternating current
- ⑦ Follow the instructions for use
- ⑧ Disposal
- ⑨ Non-ionizing radiation
- ⑩ Defibrillation-protected application part
- ⑪ Date of construction


Unit Front

- 1. Controller
- 2. Output HF CUT
- 3. Neutral electrode
- 4. Switch cut or cut/coag
- 5. LED Cut
- 6. LED Cut/COAG

Unit Back

- 1. ON/OFF switch
- 2. 2230 V Power supply with fuses
- 3. Footswitch
- 4. Label

Standard hand piece

HF-electrode CUT






Neutral electrode

Precision of power setting radio frequency (HF)

The power of the radio frequency generator depends on the tissue specific resistance and can vary within limits. The specified 50 W are in accordance with a specific resistance of 1 kΩ. The settings of the radiofrequency will be displayed without the measuring unit [W] and will be scaled in accordance with the graph.

Selection of the correct radio frequency (HF) mode

ATTENTION

All surgical applications with radio frequency need local or block anaesthesia. Make sure the tissue is moist.

Cutting - CUT MODE

This mode – with a permanent power flow – is best rated for clean cuts in tissue without coagulation. In this mode cutting is with marginal heat and little hemostasis and can be used near to bones or to the periosteal tissue to avoid shrinking processes of the tissue. Histological examinations can be done in this mode as well.

Hint: Activate the electrode by pressing the foot switch before touching the tissue. Then a regular cut will be produced from the beginning.

Cutting / Coagulation - CUT / COAG MODE

This mode allows the precise cutting and simultaneous coagulation of the cut surface. Clinically the coagulation zone is marginal, but allows an effective hemostasis, does not disturb the primarily wound healing and disappears spontaneous after the wound is healed. These cuts do not need suturing, thus this mode is very effective in cosmetic surgery.

Hint: Activate the electrode by pressing the foot switch before touching the tissue. Then a regular cut will be produced from the beginning.

Neutral electrode = NE

When operating the HF mode, it must always be worked with the connected NE. This ensures an optimal power during use. The NE must be placed between the back of the patient and the treatment chair, as close as possible to the head..

Practical exercise on a beef model

Prepare the device for operation and follow the steps mentioned below.

1. Select a piece of fresh, lean beef. Veal is not suitable, because it does not change color when cut with an electrode. Because of its cell structure pork is also not suited. Wait until the beef has obtained room temperature.

Make sure that the meat is placed on the plugged-in neutral electrode. If this is not the case the waves cannot derive and, consequently you cannot proceed with the exercise.

2. Insert the electrode of your choice (Multi Tip, loop, diamond etc.) into the hand piece.
3. Turn the intensity control dial to 8.
4. Push the toggle switch to **Cut**
5. Activate the foot switch
6. Make several incisions of different lengths and depths with even, brushing movements. Then take the energy from the electrode and look at the result. You will notice that the intensity adjustment was too high, which caused sparks and remarkable discoloration along the cutting line.
7. Reduce the intensity setting to 1. You will notice that the electrode will either cut only if it is dragged and pulled through the meat, or does not cut at all. If a cut is possible at all, bits and pieces of meat will get caught on the electrode.
8. Repeat the above described process with slightly increased settings, until you reach the point at which discolorations and visible discharges of sparks no longer appear. The tip of the electrode should not encounter resistance. The cut should be precisely even, without the occurrence of a discharge of sparks and without the necessity to drag the electrode. Continue with your endeavors by executing slow, intermediate and faster cuts for each particular setting in order to achieve the necessary expertise and confidence you will need for the surgery of your patients.
9. Turn the toggle switch to cut with simultaneous coagulation and repeat the exercises. You will see that cuts with the slightly modified wave require a higher setting than when the fully filtered wave is used. This is normal and should be taken into consideration for the later work on the patient.

Application examples: HF

HF CUT & CUT COAG				
PRG No.	Program	Power (Watt)	Coagulation grade	Indications & remarks
1	CUT	ca. 28	6	(filtered wave) - sulcus dilatation - gingivectomy - internal gingivectomy - open curetage
3	CUT	ca. 22	5	- tumor resection - lap preparation - vestibulum plastic - excision
2	CUT COAG	ca. 22	6	(slightly modulated or non filtered wave) - gingivoplastic - exposure of teeth, stubs, approximal steps or crown edges - removal of hyperplasia
4	CUT COAG	ca. 17	7	for ablation of tissue if simultaneous coagulation is requested with the cut (Attention: 10% loss of tissue about 24 h post operatively due to extended lateral heat) Use only if distance to bone or periost is sufficient!

Processing instructions according to EN ISO 17664:2021 of HF handpiece 452 423 (yellow), 452 425 (blue)

Warnings: Attention! The product has to be cleaned, disinfected and sterilized according to these processing instructions before its first use as well as after every further use. Disinfection alone is not enough. Thorough cleaning and disinfection are decisive prerequisites for the effective sterilization. The processing should be started as early as possible, 2 hours after use, however, at the latest. During the processing, the product should not be unnecessarily exposed to wetness or humidity. Please also observe the legal regulations applicable in your country as well as the hygiene regulations of the medical practice and/or the hospital.

Limitations in the processing: The frequent processing has only minor effects on this product. The end of the product's useful life is mainly determined by wear and damage caused by use. In case of doubt, the products should always be sorted out and replaced early. The user alone is responsible for making the decision regarding repeated use of the product. If the product is used too frequently, the manufacturer does not accept any warranty for the function, performance and safety of the product.

Personal protective equipment: For reasons of occupational safety and to minimize infections / cross infections, suitable personal protective equipment (protective clothes, protective gloves, protective goggles and nose/mouth protection) has to be used during the entire processing.

INSTRUCTIONS:

INITIAL TREATMENT AT THE PLACE OF USE

Remove surface contamination using a disposable cloth / paper cloth when you are still at the place of use. Rinse the products with water (at least drinking water quality) 2 hours after the application at the latest. The drying of residues or contamination of any kind on the product has to be avoided. When using the product for the first time at the place of use, do not use any aldehyde-containing or alcohol-containing agents as they may cause protein fixation.

PREPARATION BEFORE CLEANING

Disassemble the product into its individual components. Visual inspection for damage and wear.

CLEANING: MANUAL

The manual cleaning should be completed in a water bath (at least drinking water quality) with the cleaning agent specified below using a brush below the water surface to achieve both, sufficient cleaning of the products / individual components avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product / individual component; it should, however, not be less than 1 minute.

Afterwards, rinsing of the products / individual components with water (at least drinking water quality) for at least 20 seconds.

Description of the recommended cleaning process:

Pre-cleaning: Remove surface contamination using a disposable cloth

Cleaning: Insert the instrument into the strainer bowl. Hang the strainer bowl into the cleaning bath with cleaning solution, remove contamination by means of cleaning brushes.

Cleaning tools: Flexbrush brush, REF 605 254, Mirabrush, REF 605 260/61, Apply-Tips, REF 605 510

Rinsing: Rinse with water (at least drinking water quality) for 20 sec.

Drying: Drying at room temperature

The manual cleaning should not exceed a temperature of 45 °C.

It has to be ensured that visual contamination has been completely removed. If after the cleaning process, there is still contamination visible, the cleaning has to be repeated.

DISINFECTION: MANUAL

For the manual disinfection, please only use approved disinfectants with tested effectiveness (CE mark, VAH-/DGHM-listed).

Put the products / individual components into corresponding disinfection baths according to the disinfectant manufacturer's specifications. Ensure that the products are sufficiently covered and do not contact each other. The instructions for use of the disinfectant manufacturer must be strictly observed. Particularly the concentrations and contact times have to be observed.

Application with the following agents is recommended: - Dürr Dental ID 213 instrument disinfection

Description of the recommended disinfection process:

Completion of proper cleaning preparation and manual cleaning

Disinfection: Insert the part of the instrument with patient contact into the strainer bowl. Hang the strainer bowl into the disinfection bath with disinfection solution

Disinfectant: ID 213, company Dürr Dental, concentration: 2 %

Contact time: 5 min

Rinsing: Rinse with water (at least drinking water quality) for 20 sec.

Drying: Drying at room temperature

Make sure that the disinfectant is compatible with the products and cleaning agents that might be used. The pH of the disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents

(e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic / halogen hydrocarbons.

Please observe thermal restrictions that might apply considering the disinfectant manufacturer's specifications.

Afterwards, rinsing of the products / individual components with water (at least drinking water quality) for at least 20 seconds.

Before further processing, the product must be free from any residues and dry.

DRYING

If you process the products manually, the products / individual components may be dried by blowing them off with filtered, oil-free compressed air according to DIN ISO 8573-1 (medical cleanliness class) or at room temperature.

MAINTENANCE, CONTROL AND CHECK

The product does not make special requirements on maintenance. Re-assemble disassembled products / individual components. You must always carry out a visual inspection for contamination, damage, wear and deformation before and after the individual work steps. Damaged or corroded products must no longer be used. If the product / individual component is not visibly clean, the entire processing has to be repeated or the product / individual component has to be properly disposed of.

PACKAGING

The product is to be packed in suitable and standardized transparent sterile packaging (sterilization bags) and sealed. Observe the instructions of the sterilization bag and sealing machine manufacturers as well as the current, normative requirements. Products / individual parts that have not been sterilized in bags must be used immediately.

STERILIZATION

Only tested steam vacuum autoclaves may be used. Make sure that during the sterilization of several different products / individual components, the autoclave is not overfilled and that the products / individual components do not contact each other.

The following sterilization cycles may be completed:

Steam sterilization, 134 °C, hold time 5 minutes or steam sterilization, 121 °C, hold time 15 minutes

To dry the products / individual components, the drying cycle of the autoclave should be set. Use the autoclave manufacturer's instructions for use.

STORAGE

To maintain the sterility, the products have to be stored in standardized sterilization bags at a dry, clean place until they are used. If the sterile packaging is damaged, the products must be processed once again before their use.

ADDITIONAL INFORMATION

Mark the sterilized products according to the legal and national regulations. The recommended storage duration for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.

The user must ensure within the scope of their quality management system that specified processing cycles that might apply (see Limitations in the processing) are not exceeded.

The manufacturer and the competent authority of your member state shall be immediately notified all serious incidents occurring in connection with the product.

For the disposal of the products, there are regionally differing recommendations and regulations. Ask your competent disposal company for the current regulations in your region.

Instructions for Use of the RF handpieces

Inserting the electrode into the hand piece

- Turn the head of the hand piece one turn to the left.
- Insert the shaft of the electrode all the way into the hand piece.
- Turn the head of the hand piece to the right.

Now the electrode should fit tightly in the hand piece!

- Check the tightness of the electrode by gently pulling it. It has to be impossible to turn the electrode or to drag it out of the shaft.
- Make sure that the non-coated part of the electrode (metal part) is inserted completely in the hand piece. No metal parts of the shaft should be visible outside of the head of the hand piece. Otherwise, you pose a risk of injury to you and your patients.

IMPORTANT: Only the approved original electrodes of Hager & Werken must be used. These electrodes are adapted to the hand pieces and devices.

Inserting the cables of the hand pieces into the devices

The hand pieces are equipped with yellow or blue cables and matching plugs.

Put the plug into the corresponding socket, showing the same colour the cable and cable plug have.

Pay attention to inserting the plug very deep – all the way to the stop.

The following applies to the hf1 Surg bipolar only

Always put both hand pieces into the colour-coded sockets and according receptacles. This is the only way to use the device. By removing one hand piece from the receptacle, the second hand piece becomes functionless.

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Environmental Protection Directive

According to the EU Directive 2012/19/EU, the following applies to the disposal of electronic and electrical equipment:

These devices must not be disposed of with household waste .

At the end of its service life, the user is legally obliged to take this equipment to the public collection points set up for this purpose or to return it to the collection points set up for this purpose or to return them to the point of sale.

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Copyright © 2024 - Hager & Werken GmbH & Co KG
HAGER & WERKEN GmbH & Co. KG
Arazy Group Swiss GmbH