

POLISHING INSTRUMENTS
FOR THE DENTAL PRACTICE

ROTA
DIATECH





1 COMPREPOL ULTRA				
ISO Ø 1/10 mm	030	050	060	100
Head L in mm	7.5	10	7.5	2.5
Polisher type	2101RA	2103RA	2111RA	2127RA
RA 10 pcs	230055AA	230059AA	230063AA	230069AA
RA 50 pcs		230060AA	230064AA	

2 COMPOSHINE ULTRA				
	030	050	060	100
	7.5	10	7.5	2.5
	2201RA	2203RA	2211RA	2227RA
	230072AA	230076AA	230080AA	230086AA
		230077AA	230081AA	

Composite, Compomere



1 COMPREPOL PLUS						
ISO Ø 1/10 mm	030	040	060	100	110	140
Head L in mm	7.5	10	7.5	2.5	1.6	1.6
Polisher type	2301RA	2303RA	2311RA	2327RA	235G11RA	235G14RA
RA 5 pcs	230291AA	230293AA	230295AA	230297AA	60019877	60019878
RA 10 pcs		230294AA	230296AA			

2 COMPOSHINE PLUS						
	030	040	060	100	110	140
	7.5	10	7.5	2.5	1.6	1.6
	2401RA	2403RA	2411RA	2427RA	245G11RA	245G14RA
	230299AA	230301AA	230303AA	230305AA	60019879	60019880
		230302AA	230304AA			

Composite, Compomere



1 PREBROWN			
ISO Ø 1/10 mm	040	045	060
Head L in mm	7.5	10	9.5
Polisher type	1101RA	1103RA	1111RA
RA 10 pcs	230001AA	230005AA	230013AA
RA 50 pcs		230006AA	

2 HIGREEN			
	040	045	060
	7.5	10	9.5
	1201RA	1203RA	1211RA
	230028AA	230032AA	230040AA

DIASHINE COMPOMANT PLUS				
	030	040	060	100
	7.5	10	7.5	2.5
	9101RA	9103RA	9111RA	9121RA
5 pcs	230264AA	230265AA	230267AA	230271AA
10 pcs		230416AA		

Precious Metal Alloy, Amalgam

Composite, Compomere




1 ZIRCOPOL PLUS			
ISO Ø 1/10 mm	040	060	140
Head L in mm	10	7.5	1.6
Polisher type	9703RA	9711RA	975G14RA
RA 5 pcs	230379AA	230381AA	60022004
RA 10 pcs			

2 ZIRCOSHINE PLUS			
	040	060	140
	10	7.5	1.6
	9803RA	9811RA	985G14RA
	230391AA	230393AA	60022005

UNISHINE			
	055	070	090
	15	10	8.5
	4203RA	4211RA	4216RA
	230203AA	230211AA	230215AA

CAD/CAM Restoration, Zirconia

Universal



1 CERAPRESHINE		
ISO Ø 1/10 mm	050	070
Head L in mm	10	10
Polisher type	3201RA	3211RA
RA 5 pcs		
RA 10 pcs	230114AA	230122AA

2 CERASHINE		
ISO Ø 1/10 mm	050	070
Head L in mm	10	10
Polisher type	3301RA	3311RA
RA 5 pcs		
RA 10 pcs	230139AA	230147AA

1 CERAFIN PLUS				
ISO Ø 1/10 mm	040	060	110	140
Head L in mm	10	7.5	2	1.6
Polisher type	9303RA	9311RA	9331RA	935G14RA
RA 5 pcs	230317AA	230319AA	230325AA	60019881

Ceramic

Ceramic



2 CERAPRESHINE PLUS					
ISO Ø 1/10 mm	030	040	060	110	140
Head L in mm	7.5	10	7.5	2	1.6
Polisher type	9401RA	9403RA	9411RA	9431RA	945G14RA
RA 5 pcs	230331AA	230335AA	230337AA	230343AA	60019882

3 CERASHINE PLUS					
ISO Ø 1/10 mm	030	040	060	110	140
Head L in mm	7.5	10	7.5	2	1.6
Polisher type	9501RA	9503RA	9511RA	9531RA	955G14RA
RA 5 pcs	230349AA	230353AA	230355AA	230361AA	60019883

Ceramic

Ceramic



BRUSH			
ISO Ø 1/10 mm	060	040	040
Head L in mm	5.5	6.5	6.5
Polisher type	7001RA	7002RA	7003RA
RA 5 pcs	232002AA	232004AA	232006AA
RA 10 pcs	232001AA	232003AA	232005AA

MANDREL	SWISSFLEX DISCS				
Length	Grit size	C/70µm	M/50µm	F/30µm	UF/5µm
25.4	Pack size	80 pcs	80 pcs	80 pcs	80 pcs
702RA	Polisher type	704UM ¹ /703UM ²	706UM ¹ /705UM ²	708UM ¹ /707UM ²	710UM ¹ /709UM ²
232018AA	Head Ø: 9 mm	232020AA ¹	232022AA ¹	232024AA ¹	232026AA ¹
	Head Ø: 13 mm	232019AA ²	232021AA ²	232023AA ²	232025AA ²

Ceramic, Composite, Metal

Composite, Amalgam, Glass Ionomer Cement, Semi-precious Metal, Precious Metal

SWISSFLEX STRIPS	
Length × Width	175 mm × 3.5 mm
Pack size	80 pcs
700UM: C(70µm)/M(50µm)	232016AA
701UM: F(30µm)/UF(5µm)	232017AA

See SwissFlex Discs

DIAStrip	
Length × Width	80 mm × 2.5 mm
Polisher type	600UM
5 pcs	232014AA
10 pcs	232015AA

Composite, Compomere, Glass Ionomer Cement, Amalgam



Order information

ISO number (material, shank, shape, grit size, variable diameter) ————


Head diameter in 1/10 mm ————

Head length ————

Polisher description ————

Packaging sizes available ————

ISO 658 204 243 523 ...



Illustration

ISO Ø 1/10 mm	050
Head L in mm	10
Polisher type	2103RA
RA 10 pcs	230059AA
RA 50 pcs	230060AA

REF in respective packaging size

Head Shapes

ISO Head Shape	Description	ISO Head Shape	Description
014	Small Wheel	243 040 10	Point
030	Cup	243 045-055 010	Flame
042	Wheel	273	Dome
044	Wheel	274	Bullet
107	Cylinder	292	Torpedo
114	Cylinder	303	Large Lens
173	Taper	304	Small Lens
243 030-040 7.5	Minipoint	372	Disc
		543	ShapeGuard

Recommended Speed Ranges

Polishing System	RPM	Polishing System	RPM
Comprepol Ultra	5'000–10'000	Cerapreshine	5'000–10'000
Composhine Ultra	5'000–10'000	Cerashine	5'000–10'000
Comprepol Plus	3'000–8'000	Cerafin Plus	5'000–10'000
Composhine Plus	3'000–8'000	Cerapreshine Plus	5'000–10'000
Prebrown	5'000–10'000	Cerashine Plus	5'000–10'000
Higreen	5'000–10'000	Brush	5'000–10'000
Diashine Compomant Plus	7'000–10'000	SwissFlex	10'000–16'000
	Prepolishing		
	3'000–8'000	Specialty Head Shape	RPM
	High lustre	ShapeGuard	10'000–12'000
Zircopol Plus	5'000–10'000		
Zircoshine Plus	5'000–10'000		
Unishine	5'000–10'000		

© 2017 Coltène/Whaledent AG – www.coltene.com



CH **Coltène/Whaledent AG**
Feldwiesenstrasse 20
9450 Altstätten / Switzerland
T + 41 71 757 53 00
F + 41 71 757 53 01
info.ch@coltene.com

UK **Coltène/Whaledent Ltd.**
The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF / UK
T +44 1444 235486
F +44 1444 870640
info.uk@coltene.com



COLTENE