

ProTaper Ultimate™ Solution

A solution combining the latest generation of ProTaper files, enhanced disinfection and dedicated obturation that works seamlessly together.

Shaping

STANDARD

Slider 016.002v



Shaper 020.004v



F1 020.007v



F2 025.008v



F3 030.009v



LARGER AND STRAIGHTER CANALS ONLY

Slider 016.002v



FX 035.012v



FXL 050.010v



Cleaning



SmartLite® Pro EndoActivator®

Coming soon

Irrigation Needle



Obturation


Absorbent Points



AH Plus® Bioceramic Sealer



Conform Fit® Gutta-Percha Points



ProTaper Ultimate™

Product references

ProTaper Ultimate™ Rotary Files

SX 19 mm PULRSX19 x6



Slider 21 mm PULRSL21 x6



25 mm PULRSL25 x6



31 mm PULRSL31 x6




Shaper 21 mm PULRSH21 x6



25 mm PULRSH25 x6



31 mm PULRSH31 x6



F1 21 mm PULRF121 x6



25 mm PULRF125 x6



31 mm PULRF131 x6



F2 21 mm PULRF221 x6



25 mm PULRF225 x6



31 mm PULRF231 x6




F3 21 mm PULRF321 x6



25 mm PULRF325 x6



31 mm PULRF331 x6



FX 21 mm PULRFX21 x6



25 mm PULRFX25 x6



31 mm PULRFX31 x6



FXL 21 mm PULRFXL21 x6



25 mm PULRFXL25 x6

31 mm PULRFXL31 x6

ASSORTMENT FILE PACKS

Sequence 21 mm PULRSQC21 x5

Slider, Shaper, 25 mm PULRSQC25 x5

F1, F2, F3 31 mm PULRSQC31 x5

ProTaper Ultimate™ Manual Files

SX 19 mm PULHSX19 x6



Slider 21 mm PULHSL21 x6



25 mm PULHSL25 x6



31 mm PULHSL31 x6



Shaper 21 mm PULHSH21 x6



25 mm PULHSH25 x6



31 mm PULHSH31 x6



F1 21 mm PULHF121 x6



25 mm PULHF125 x6



31 mm PULHF131 x6



F2 21 mm PULHF221 x6



25 mm PULHF225 x6



31 mm PULHF231 x6



F3 21 mm PULHF321 x6



25 mm PULHF325 x6



31 mm PULHF331 x6



FX 21 mm PULHF21 x6



25 mm PULHF25 x6



31 mm PULHF31 x6



FXL 21 mm PULHFXL21 x6



25 mm PULHFXL25 x6

31 mm PULHFXL31 x6

ASSORTMENT FILE PACKS

Sequence 21 mm PULHSQC21 x5

Slider, Shaper, 25 mm PULHSQC25 x5

F1, F2, F3 31 mm PULHSQC31 x5

Rx only

Files Dentsply Tulsa Dental Specialties, 608 Rolling Hills Drive, Johnson City, TN 37604, USA

Gutta-Percha & Absorbent Points Dentsply Tulsa Dental Specialties, 608 Rolling Hills Drive, Johnson City, TN 37604, USA

Irrigation Needle Produits Dentaires SA, Rue des Bosquets 18, CH-1800, Vevey

AH Plus® Bioceramic Sealer Dentsply Tulsa Dental Specialties, 608 Rolling Hills Drive, Johnson City, TN 37604, USA
Manufactured for

SmartLite® Pro EndoActivator® Dentsply LLC, 38 West Clarke Avenue, Milford, DE 19963 USA

ProTaper Ultimate™

Performance unlimited

ProTaper Ultimate™ Conform Fit® Gutta-Percha Points

Box of 60 gutta-percha points

F1 PULGF1

F2 PULGF2

F3 PULGF3

FX PULGFX

FXL PULGFXL

ASSORTMENT

F1-F3 PULGAST

GP points assortments: 12x - F1, 24x - F2, 24x - F3

ProTaper Ultimate™ Absorbent Points

Box of 180 absorbent points

F1 PULPF1

F2 PULPF2

F3 PULPF3

FX PULPFX

FXL PULPFXL

ASSORTMENT

F1-F3 PULPAST

PP points assortments: 30x - F1, 90x - F2, 60x - F3

AH Plus® Bioceramic Sealer



Starter Kit AHPBIOSEAL

Syringe Refill AHPBIOREFILL

Cannula Refill AHPBIOCANNULA

dentsplysirona.com/protaperultimate
#endopassion




NAM: Technical modifications on our products are not subject to notifications. Photos of our instruments are not contractual. © 2021 Dentsply Sirona Inc. Rx Only. INL / T EN PTUL BRC 000 / 00 / 2021 - created 07/2021


Shaping

1 Slider-Shaper-Finishers sequence to cover a full range of anatomical situations

Deep Shape by 



Can be used directly to start the glide path in the majority of cases¹



Enhanced cutting efficiency and hauling of debris in the coronal two-thirds



In the apical third, give you the flexibility for the most challenging anatomies, without breaking or unwinding




Our difference

100% of users consider that ProTaper Ultimate™ files provide a predictable root canal preparation¹

Shaping time



Stronger for longer endurance



60% faster than Komet F360 and Brasseler EndoSequence® CM
130% faster than EdgeEndo® EdgeTaper Platinum™²

¹ Based on a 2021 user evaluation of dental clinicians. Data on file.

² Internal testing. Data on file.

All product and company names are trademarks™ or registered® trademarks of their respective holders. Use of them does not imply any affiliation with or endorsement by them.


Cleaning

More than 25% increased cleaning efficacy vs. sodium hypochlorite soak without activation¹


Irrigation Needle



The irrigation needle **curves and flex easily** providing an efficient cleaning and disinfection until the apex





EndoActivator




Activated fluids promote deep cleaning and disinfection facilitating 3-D obturation and long-term success



Our difference



Comparison of mean scores for lateral canals obturated



The EndoActivator® provided better obturation of lateral and accessory canals and resulted in less remaining debris*

¹ Data on file, 2021.

* A quantitative and qualitative analysis of ultrasonic versus sonic endodontic systems on canal cleanliness and obturation Valerie Kanter, DMD, Emily Weldon, DMD, Uma Nair, DMD, MDS, Claudio Varella, DDS, MS, Keith Kanter, DDS, Kenneth Anusavice, PhD, and Roberta Pileggi, DDS, MS, Gainesville, FL UNIVERSITY OF FLORIDA. (Oral Surg Oral Med Oral Pathol Oral Radiol Endod 2011;112:809-813)

Obturation

Obturation with an ideal fit 100% of the time

Bioceramic Sealer



Biocompatible and does not stimulate periodontal tissue



Conform Fit® Gutta-Percha Point




Less waste and saves time because the first cone fits everytime





Our difference

Setting time



10 hours EndoSequence® BC Sealer™



Range of fit with Conform Fit® Gutta-Percha Point

21% fit¹ When using traditional gutta-percha cones



Range of fit with Traditional gutta-percha

¹ Data on file, 2021.