

CE0120



Committed to Excellence  
**Technical & General Ltd**


Code: KF 00

## K- Files




- Exceeds all ISO/ADA specifications
- High Quality Stainless Steel
- Hand use with endo stops
- Blister card packing 6 pieces/blister in lengths of 21 , 25 , 28 & 31 mm
- Colour coded, ISO colours & sizes



- Cross section  ISO 06-25 Square X section for better resistance to fracture

- Cross section  ISO 30-140 triangular for better resistance to fracture

- Shape 

- Handlee **CC-Cord® Griff** **CC-Cord® handle** **Symbol**  
  



- Excellent cutting ability thanks to Cr-Ni Steel with high cutting ability, precision cutting edges and unique grinding technique
- Flatter cutting angles towards the tip ensure high resistance to fracture
- Special steel alloy with Uniform structure for maximum resistance to fracture & flexibility
- The inactive tip conducts the instrument safely along the canal
- Sterilizable up to 200 C /390 F

### Available Types & Sizes

Length	Size
21mm	06, 08, 10, 15, 20, 25, 30, 35, 40, 15-40, 45, 50, 55, 60, 70, 80 & 45-80
25mm	06, 08, 10, 15, 20, 25, 30, 35, 40, 15-40, 45, 50, 55, 60, 70, 80 & 45-80
28mm	15, 20, 25, 30, 35, 40, 15-40 & 45-80
31mm	15, 20, 25, 30, 35, 40, 15-40 & 45-80

### Sales and Information

Technical & General Ltd., 2 Albion Place, London W6 0QT  
 Tel. 020 8834 7575 Email: sales@technical-general.com  
 www.dentistmarket.com www.tgdent.com

