

UPDATE No. 22
January 2010

C

a

v

e

x

N

e

w

s

NEW! The Cavex VacuFormer System: From impression to the desired splint/bleaching tray/mouthguard in 1 session..!

The Cavex VacuFormer System is a complete vacuum thermo forming system comprising the Cavex VacuFormer, laminates, a unique plaster composite and accessories. The Cavex VacuFormer System enables the dentist to optimize the treatment of and service to the patient. Bleaching trays, splints, bruxism guards and mouth protectors can be prepared during patient visits without the intervention of the dental laboratory. A detailed protocol, illustrated with photographs, enables anyone to master this technique.

The Cavex VacuFormer is a “no-nonsense” all-round thermal vacuum former. This robust machine has an electric, high performance vacuum motor suitable for every conventional dental deep-drawing application. The Cavex VacuFormer is easy to operate and is protected against overheating.



Cavex Instant Stone (art. nr. VT010) is a unique plaster composite that can cast a perfect plaster mould in just 10 minutes. Thanks to a sophisticated blending system, Cavex Instant Stone is very clean and easy to use.



The Cavex Mouth Protectors (art. Nr. VT020) are professional 3.8 mm dual layered laminates with a combined hard and soft layer. The soft layer ensures good shock absorption and the hard layer distributes the force of the impact, resulting in optimal protection of your teeth.



UPDATE No. 22
January 2010

C
a
v
e
x

N
e
w
s

The Cavex Bleach laminates (art. Nr. VT022) are made of a soft, clear material that is easy to cut and finish. The soft material is comfortable to the gingiva and thus ideal for making bleaching trays.



The Cavex Splint laminates (art. nr. VT021) are made of a transparent and rigid material which is ideal for splints, surgical trays, orthodontic retainers, etc. This unique, clear material makes the Cavex Splint invisible when worn in the mouth.



The Starter Kit (art. nr. VF001) consists of all equipment and materials that are required for making ready-made deep-drawing applications:

1 Cavex VacuFormer + 18 jars
Instant Stone + laminates
(50 x Splint + 25 x Bleach +
12 x Mouth Protector) +
Scissors + 6 Base trays +
protocols per indication.

