



BETTER TOOLS FOR BETTER RESULTS

At Ultradent, we appreciate your precise craftsmanship and know that your hard work enhances patient comfort and helps to improve oral health. Our quality line of lab products is designed to make your job easier from start to finish.



ULTRADENT
PRODUCTS, INC.

800.552.5512 ULTRADENT.COM




DENTAL LAB PRODUCTS AND TOOLS

Ultradent® LC Block-Out Resin

Ultradent LC Block-Out Resin is adept at filling in voids and blocking out undercuts. With its high contrast, optimal viscosity, and syringe delivery, it is easy to place accurately and precisely. It is also ideal for creating reservoir space in whitening trays. Studies show that reservoirs can improve whitening results by keeping more gel in contact with teeth.

Provides reservoir space for whitening trays and is useful for other laboratory procedures such as model and die repairs.

- 240 **Ultradent LC Block-Out Resin Kit**
4 x 1.2 ml Ultradent LC Block-Out syringes
20 x Black Mini tips



- 242 **Ultradent LC Block-Out Resin Econo Kit**
20 x 1.2 ml Ultradent LC Block-Out syringes
20 x Black Mini tips



- 241 **Ultradent LC Block-Out Resin Refill**
4 x 1.2 ml Ultradent LC Block-Out syringes



Upon curing, the strong, no-mix material becomes hard instantly—perfect for blocking out undercuts on dies and filling in voids.

Sof-Tray® Classic Sheets

Select the 0.035" for most whitening trays, the 0.060" for whitening patients who are bruxers, and the 0.080" for severe bruxers or as a TMJ appliance.

- 226 **Sof-Tray Regular Sheets - 0.035" 5"x 5"**
25 x Tray sheets (0.9 mm - 127 x 127 mm)

- 227 **Sof-Tray Medium Sheets - 0.060" 5"x 5"**
20 x Tray sheets (1.5 mm - 127 x 127 mm)

- 284 **Sof-Tray Heavy Sheets - 0.080" 5"x 5"**
20 x Tray sheets (2.0 mm - 127 x 127 mm)



Ultradent® Utility Vinyl Cutters

Ideal for gross trimming of tray.

- 604 **Ultradent Utility Vinyl Cutters**



Ultradent® Ultra-Trim Scalloping Scissors

- For precision trimming of border around interdental papilla
- Spring loaded to minimize finger fatigue
- Grip tray material easily
- Made of durable stainless steel

- 605 **Ultradent Ultra-Trim Scalloping Scissors**

