

#whdentalwerk



wh.com



System overview adaptor

Assistina Twin, Assistina One & Q-Dry

Straight and contra-angle handpieces

All straight and contra-angle handpieces with ISO connection

Contra-angles with removable heads

T1 Classic from Sirona Dental Systems for contra-angle handpieces

FG chuck system (W&H, Sirona Dental Systems)

W&H handpiece and contra-angle head



Adaptor REF 02679000



Adaptor REF 02691000



Adaptor REF 02693000



Adaptor REF 05204600

Adaptor Quick ISO REF 06422600



Turbines

W&H turbines, air scaler and air motors with Roto Quick system

T1, T2, T3 from Sirona Dental Systems

Turbines with Multiflex connection from Kaltenbach & Voigt

Turbines with Borden 2(3)-hole connection

Turbines or adaptors with standard 4-hole connection (with or without light)



W&H Adaptor REF 02690400



Adaptor REF 02692000



Adaptor REF 04713200



Adaptor REF 02083500

Adaptor Quick RM REF 06395100

