

X-Smart® Apex Locator

The right way to the apex





X-Clusive technology for X-Treme precision

X-Smart® Apex Locator offers highly reliable accuracy—with patented measurement technology, an easy-to-use touchscreen interface, and audio/visual cues for added guidance during procedures.



6



**HIGHLY
PRECISE**
working length
determination*



Measuring with X-Pert accuracy

Other apex locators use mixed AC signal frequencies to detect the apex; X-Smart® Apex Locator alternates them, canceling the need for signal filtering and eliminating the noise that can interfere with performance.



Get up to **150**
root canal
treatments
with a full
charge*



Features designed for you

- 1 Large, full-color, high-resolution touchscreen that works with gloves
- 2 Dr's Choice feature allows for saved user preferences and provides visual and audio cues when the file tip has reached the pre-selected position
- 3 Visual and audio guidance during measurements
- 4 Long-lasting, lithium-ion rechargeable battery with fast charging
- 5 USB Type-C connector for more durability and flexible orientation
- 6 Connection of File Clip and Lip Clip to any socket of the measurement cable with onboard back storage

An accurate working length determination is crucial to the success of your root canal treatment.


Let X-Smart® Apex Locator guide your way. Learn more at dentsplysirona.com/XSmartApexLocator

Technical Specifications

- 5" touchscreen interface
- Lithium-ion battery
- Fully recharged in approximately 3 hrs
- On/Off button
- Sound control
- Device settings
- Automatic turn-off after 5-15 min
- Dimensions: Width 95 mm x Depth 146 mm x Height 111 mm
- Calibration not required
- Frequencies used for measurements 500 Hz and 8 kHz
- Product Warranty: 24 months / Accessories Warranty: 6 months



Product Reference

X-Smart® Apex Locator		ACCESSORIES		
	B00XSAL000000	1x X-Smart® Apex Locator	B00XSALAC1000	X-Smart® Apex Locator Measurement Cable
		1x Measurement Cable	B00XSALAC2000	X-Smart® Apex Locator File Clip, set of 2
		2x File Clip	B00XSALAC3000	X-Smart® Apex Locator Lip Clip, set of 5
		2x Lip Clip	B00XSALAC4000	X-Smart® Apex Locator Touch Probe, set of 2
		1x Touch Probe	B00XSALAC5000	X-Smart® Apex Locator Charger with Plugs (EUR, UK, US, AUS)
		1x Charger with Plugs (EUR, UK, US, AUS)		